


HKS880B CAVITATION+RADIO FREQUENCY



2 Probe 25K+40K

APPLICATION

- *Intensive physical lipolysis to remove surplus fat
- *Burn surplus fat, Lymphatic Treatment, skin tightening, improve skin elasticity
- *Improve orange peel organization, body shaping


PARAMETER

- *Energy:3-6W/cm2 can adjust
- *Voltage:110V or 220V
- *Power:150W
- *Frequency:25K+40KHz
- *Active Surface:20cm2
- *Operate timer:1-20min
- *Display:8.0inch touch screen
- *Emission model:Scan 2% compared with autotuning Frequency
- *Functionality pressure:500mbar-1060mbar
- *Certifications:CE LVD,CE EMC
- *Package size:42*37*33cm
- *N.W:5.5kg
- *G.W:9.0kg

2 Probe 3MHZ

APPLICATION

- *Skin lifting
- *Improve conditions of sagging skin
- *Forehead wrinkle removal
- *Breast lifting
- *Shape body contour
- *Hand care, effective on removal wheelk and scars



Inter face

